# atemious<sup>®</sup> PRO2XXS

MADE IN GERMANY

FFP2 VIREN-FILTERNDE ATEMSCHUTZMASKE  
VIRUS FILTERING RESPIRATOR MASK

ATEMIIOUS PRO2XXS  
ART. 2004  
EN 149:2001+A1:2009  
FFP2 NR

CE 2233

UNIVENT<sup>+</sup>  
MEDICAL

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*We thank you for your trust, with the **atemious pro2xxs** you have chosen a German quality product. Univent GmbH has been in existence for more than 25 years and manufactures high quality respiratory protection masks. The products are 100% made in Germany and meet the highest quality standards according to CE guidelines. Univent Medical has its own test laboratory, where the superior filter quality is continuously controlled and ensured.*

Univent Medical GmbH  
Auf Luckenburg 12  
78056 Villingen-Schwenningen  
[univent-medical.de](http://univent-medical.de)

**UNIVENT**<sup>+</sup>  
MEDICAL

The foldable FFP2 respirator **atemious pro2xxs** offers high wearing comfort and effective protection of the respiratory tract against solid and liquid aerosols according to EN 149: 2001 + A1: 2009.

### INSTRUCTIONS FOR USE

Before using the mask, the following information must be read and followed. Failure to follow the instructions for use can result in illness or death!

- The user of the mask must be instructed on how to use it correctly.
  - Masks do not protect against gases and vapors.
  - The oxygen content in the breathing air must be minimum 19.5% by volume.
  - Respiratory protection must not be used if the concentration, type or characteristics of the harmful substances are not known.
  - The mask must be changed immediately if it is damaged or the breathing resistance increases.
- After the end of a working shift, the mask must always be changed and is not reusable (only to be used for one shift)
- Never make changes or modifications to the mask.

### INSTRUCTIONS FOR PUTTING ON THE RESPIRATOR

Wash hands before putting on the mask. Prevent contact with the inside of the mask.

Please note:  
Only stretch the rubber band facewards. Otherwise the ear band fixation may be destroyed.

2 Take the respirator out of the packaging, unfold it and pull the earband firmly in both directions. Inspect the mask visually for damage.

3 Place the respirator in your hand with the nose area facing up. Now direct the respirator to your face - first insert your chin, then pull the mask over your nose. Fix the side ear bands.

4 For a firm fit, it is necessary to adjust the clip in the nose area to the contour of the face.

5 Checking the tightness - cover the mask with both hands and exhale vigorously. With slight pressure build-up in the mask (exhale), no air should escape at the edges. When used correctly, the filter material moves slightly when inhaling. A tight fit must be ensured. If necessary, the straps can be shortened by knotting them.

**Warning:** Hair like beards can affect the required firm fit of the mask and can cause leaks.

### STORAGE AND SHELF-LIFE

Shelf life when properly stored max. 3 years. Check expiry date on masks before use. Storage and transport clean and dry between +10°C to +25°C at max. 60% humidity in original packaging. Protect against UV radiation!

### DISPOSAL

In the event of damage, high pollution or if the breathing resistance is significantly increased, the mask should be changed and disposed in accordance with regional regulations.

### CONFORMITY ASSESSMENT

The particle filtering half mask meets the requirements of the European Directive (EU) 2016/425. Type examination and conformity assessment procedures were carried out with the participation of the following notified body: GÉPTE SZT Termelészközök Felülvizsgáló és Karbantartó Kft. . Jablonka u.79. . 1037 Budapest, Ungarn

### CE 2233

Download the declaration of conformity from [www.univent-medical.de/atemiouspro2xxs/conformity](http://www.univent-medical.de/atemiouspro2xxs/conformity)

### INFO

We are happy to help you choosing the right respirator. The selection of the most suitable personal protective equipment (PPE) depends on the respective situation and may only be carried out by competent persons. For questions and further information please contact Univent Medical.

[www.univent-medical.de/info](http://www.univent-medical.de/info)